

Weekly Eyewash Activation Log

How to use this log:

- If you have more than one eyewash, use a different sheet for each location and add location descriptions as the file title or in the location field for print-outs.
- Add your initials to document completion of eyewash activation.
- Store this log with your digital Lab Safety Plan so that anyone in your lab could access it or print it out and place by eyewash station.

Learn more about emergency eyewashes and showers: [DRS Library](#) / [DRS Video](#)

Eyewashes must be activated weekly for a period long enough to verify operation and ensure that flushing fluid is clear. This helps clean out any rust, scale deposits, or bacteria that may accumulate.

Safety note: Ensure that the water will end up in a nearby drain or use a capture pan to collect water.

Eyewash Location:

Date	Initials	Date	Initials	Date	Initials	Date	Initials